## Dee Fog®

## **Anti-Fog Treatment**

Dee Fog Anti-Fog Treatment provides lens defogging while retaining clear, undistorted vision. Dee Fog's formula is compatible with mirrors or lenses, including front surfaces. It is available in liquid and spray form depending on the needs of the practitioner. Dee Fog Liquid is applied by dripping a small amount of liquid and spreading it evenly over entire mirror/lens surface. Dee Fog Spray is sprayed directly to the surface. Dee Fog Liquid is sold as a 1 oz bottle. Dee Fog Spray is sold as a 2 oz manual pump spray.

Description	Size	Item No.
Dee Fog® Anti-Fog Treatment - Spray	2 oz	0504
Dee Fog® Anti-Fog Treatment - Liquid	1 oz	0501



## **Key Features**

- No sponge needed
- Provides clear, undistorted vision
- Pleasant aroma and taste
- Convenient drip or spray application
- · No pouring and dipping required

- · Works in presence of water
- Direct application eliminates cross contamination
- No need for multiple applications
- Concentrated formula
- Manual pump spray

## **Benefits**

- Allows user to see clearly without distortion
- Eliminates cross-contamination

• Won't damage mirrors, lens, etc.

Cetylite Inc. 9051 River Road Pennsauken, NJ 08110 856.665.6111 marketing@cetylite.com

©2015 Cetylite Industries, Inc. All rights reserved. Information is summary in nature and subject to change. Dee Fog and Cetylite are registered trademarks of Cetylite Industries, Inc.

