

INTRODUCING

AccuDent[®] XD

Alginate Impression System



A beautiful smile makes
a great impression...

...but it takes a great impression to make a beautiful smile.

*Introducing AccuDent XD, the next generation
Alginate Impression System.*

AccuDent® XD


Alginate Impression System

The secret behind the system rests in the unique relationship between the two impression materials. The irreversible hydrocolloid tray and syringe gels have different viscosities, but are chemically compatible and are designed to work harmoniously in the mouth to create impressions with outstanding accuracy and detail.

The AccuDent XD Impression System creates impressions with dimensional accuracy and freedom from air bubbles and voids, making it possible to create visibly better primary stone or virtual models. These impressions are highly desirable for the fabrication of custom impression trays, cast metal frameworks, orthodontic appliances, opposing models, removable dentures, and many other types of dental cases.

Introducing AccuDent XD, the next generation Alginate Impression System.

- **Dual-phase alginate materials** for more accuracy and precision than traditional alginates
- **Thermoplastic heat moldable impression trays** to accommodate all arch forms
- **Low cost and fast set time** for optimum productivity



Additional benefits

- Color Changing properties for precise timing of impression taking
- Disposable trays to prevent cross contamination and to meet high OSHA standards
- Includes Tray Adhesive for optimal bond between tray and alginate



Delivery Forms / Product Codes

AccuDent XD Syringe Material - 12pk	673470	AccuDent XD 12 pk Edentulous 10 – (Up, Md)	673713	AccuDent XD Accessory Kit:	679893
AccuDent XD Tray Material - 12pk	673471	AccuDent XD 12 pk Edentulous 11 - (LR, Sm)	673714	(2 Syringes, 2 Soft Tips, Water Measuring Vial, 2 Finger Grip Extenders)	
AccuDent XD 12 pk Dentate 1 - (LR, Lg)	673704	AccuDent XD 12 pk Edentulous 12 - (Up, Sm)	673715		
AccuDent XD 12 pk Dentate 2 - (Up, Lg)	673705	AccuDent XD 24 pk Dentate Trays - (4 x 3ea.)	673716		
AccuDent XD 12 pk Dentate 3 - (LR, Md)	673706	AccuDent XD 24 pk Edentulous Trays - (4 x 3ea.)	673717	AccuDent XD Edentulous Trial Kit:	680375
AccuDent XD 12 pk Dentate 4 - (Up, Md)	673707	AccuDent XD Tray Adhesive - (30 ml)	673469	(24 pk Edentulous Trays, 1 ea. - Syringe Material, Tray Material, Tray Adhesive, Accessory Kit)	
AccuDent XD 12 pk Dentate 5 - (LR, Sm)	673708	AccuDent XD Syringe	679890		
AccuDent XD 12 pk Dentate 6 - (Up, Sm)	673709	AccuDent XD Syringe Soft Tips - (2 pc)	679891		
AccuDent XD 12 pk Edentulous 7 - (LR, Lg)	673710	AccuDent XD Syringe Med Tips - (2 pc)	679892	AccuDent XD Dentate Trial Kit:	680376
AccuDent XD 12 pk Edentulous 8 - (Up, Lg)	673711	AccuDent XD Syringe Finger Grip Extender	679889	(24 pk Dentate Trays, 1 ea. - Syringe Material, Tray Material, Tray Adhesive, Accessory Kit)	
AccuDent XD 12 pk Edentulous 9 - (LR, Md)	673712	AccuDent XD Water Measuring Vial	679888		

100% CUSTOMER SATISFACTION
GUARANTEED!



ivoclarvivadent.com

For more information, call us at 1-800-533-6825 in the U.S., 1-800-263-8182 in Canada.
©2015 Ivoclar Vivadent, Inc. Ivoclar Vivadent and AccuDent are trademarks of Ivoclar Vivadent, Inc.
681984 Rev. 1 8/2015

ivoclar
vivadent
passion vision innovation