

NEW

AntiSplatr® FLEX Slim

Disposable Prophy Angle

Smart Flex Technology Adjusts to Every Tooth Shape

Patent Pending



90°

Ergo






AntiSplatr® FLEX Slim

Disposable Prophy Angle




Flexing Cup

Easily adapts to the contours of any tooth




Splatter Fender

Effectively removes accumulated paste and bodily fluids



Optimal Cup Lip

Allows greater flare & subgingival access



Spiral Blades

Allow for optimum paste retention



42% Slimmer Head


Provides better visibility and access



Available in 90° and Ergo versions!

Ordering Information

- 212101 AntiSplatr Flex Slim, 200/box
- 212102 AntiSplatr Flex Slim, bulk pack, 200/box, 6 boxes/case
- 212103 AntiSplatr Flex Ergo Slim, 200/box
- 212104 AntiSplatr Flex Ergo Slim, bulk pack, 200/box, 6 boxes/case



AntiSplatr Workflow



AntiSplatr® Flex Slim
Disposable Prophy Angle

ProMate™ CL
Cordless Hygiene Handpiece

ProPaste™ One
Paste & Ring Combo





FOLLOW US



