



Thinking ahead. Focused on life.

New! TORQTECH Ultra mini

1:5 High-Speed Contra Angle

TORQTECH attachments are built tough and durable, yet provide a steady cut for even the most delicate of margins. The latest addition to the product line, the new Ultra mini, improves access and visibility with a compact head, shorter than the original, as well as competitive units. This attachment is ideal for patients with limited jaw mobility, children, and narrow treatment areas. Enjoy the stability you've come to trust, now in a smaller profile!

- Mini contra angle with 1:5 speed increasing gear ratio & optics
- Accepts standard & short shank burs
- Compact, shorter head height offers excellent access & visibility
- Low cost per use due to high cutting capacity & durability
- Newly designed water spray provides consistent cooling with four points of chip air
- Compatible with air & E-type electric motors



TORQTECH 1:5 Attachment

NEW! TORQTECH 1:5 Ultra mini Attachment



Ordering Information

TORQTECH Ultra mini: 16-1908231