



WE CARE LIKE FAMILY™

IV ADMIN & EXTENSION SETS



- Available in a variety of lengths and configurations for optimal adaptability
- For use with gravity administration sets only; not compatible with infusion pumps
- Sterile and individually packaged for optimal hygiene
- Non-pyrogenic; Not made with natural rubber latex or DEHP

ADMINISTRATION SETS

- Designed to help deliver fluids and medications safely and effectively
- Connectors can sustain a pressure of 30 seconds under 7 kg/sq. cm
- 50/bx



EXTENSION SETS

- Constructed from premium plastic materials for improved reliability
- Designed to add length to IV delivery systems
- Tubing can sustain static tensile force of 15 N for 15 seconds
- 100/bx



IV ADMINISTRATION SETS

10 drops 15 drops 20 drops 60 drops

Item	Drops/mL	Length	Needleless Port	Pre-Pierced Port	Vented	Micron Filter	Regulator	Notes	Illustration
7029	10	98" (248 cm)	2		✓		✓		
7030	10	83" (211 cm)	1	1					
7031	15	83" (211 cm)	1	1					
7032	15	78" (198 cm)		1					
7033	15	106" (269 cm)		2					
7034	15	73" (185 cm)		1					
7037	15	92" (234 cm)		2					
7041	10	89" (226 cm)		2				Primary Set with Secondary/Piggyback Port and Check Valve	
7042	10	112" (284 cm)		2				Primary Set with Secondary/Piggyback Port and Check Valve	
7044	20	92" (234 cm)		1	✓				
7045	20	83" (211 cm)		1	✓	.02	✓		
7049	20	92" (234 cm)		1	✓	.02	✓		
7051	60	83" (211 cm)		1					
7052	60	78" (198 cm)		2	✓				
7068	20	106" (269 cm)						Compatible with Zyno Pump	

IV EXTENSION SETS (NEEDLE-FREE LUER LOCK CONNECTOR)

Item	Length
7060	6" (15.24 cm)

Item	Length
7061	7" (17.78 cm)

Item	Length
7062	8" (20.32 cm)



WE CARE LIKE FAMILY™